

For Research Use Only. Not For Use In Diagnostic Procedures

Version 2.0

糖原检测试剂盒

Glycogen Assay Kit



Do not eat Store at -20° C & in the dark.

Cat.No.MKG3985

Size : 100 tests

Technical literature is available at : www.mesgenbio.com

E-mail MesGen Technical Services if you have questions on use of this system : tech@mesgenbio.com

Application

For quantitative determination of glycogen and evaluation of drug effects on glycogen metabolism.

Description

Glycogen is a branched polysaccharide of glucose units linked by α -1,4 glycosidic bonds and α -1,6 glycosidic bonds. It is stored primarily in the liver and muscle, and forms an energy reserve that can be quickly mobilized to meet a sudden need for glucose. The most common glycogen metabolism disorder is found in diabetes, in which, due to abnormal amounts of insulin, liver glycogen can be abnormally accumulated or depleted. Genetic glycogen storage diseases have been associated with various inborn errors of metabolism caused by deficiencies of enzymes necessary for glycogen synthesis or breakdown. Simple, direct and automation-ready procedures for measuring glycogen concentrations find wide applications in research and drug discovery. MesGen Biotechnology ' glycogen assay uses a single Working Reagent that combines the enzymatic break down of glycogen and the detection of glucose. The color intensity of the reaction product at 570nm or fluorescence intensity at $\lambda_{ex}/\lambda_{em} = 530/585$ nm is directly proportional to the glycogen concentration in the sample. This simple convenient assay is carried out at room temperature and takes only 30 min.

Key Features

Sensitive and accurate. Use as little as 10 μ L samples. Linear detection range in 96-well plate: 4 to 2000 μ g/mL glycogen for colorimetric assays and 0.4 to 200 μ g/mL for fluorimetric assays. Simple and high-throughput. The procedure involves addition of a single working reagent and incubation for 5-10 min at room temperature.

Samples

Serum, plasma, urine, saliva, milk, culture medium, food, agriculture etc

Kit contents

Reagent A : 5 mL Reagent B : 5 mL

Reagent C : 2 mL Reagent D : 200 μ L

20mg/mL Glycogen Standard: 1 mL

Procedures

Sample Preparation:

Samples can be prepared according as follow.

Briefly, homogenize tissue/cell sample in 25 mM citrate, pH 4.2, 2.5 g/L NaF on ice. Centrifuge 14,000 g for 5 min to remove debris, and use 10 μ L clear supernatant for the assay.

Colorimetric Procedure:

1. Equilibrate all components to room temperature. During experiment, keep thawed reagents in a refrigerator or on ice.
2. Standards and samples: Dilute standard by mixing 30 μ L Standard with 270 μ L dH₂O to give 2000 μ g/mL standard.

Dilute standard in dH₂O as follows.

NO.	2000 μ g/mL STD + H ₂ O	Vol (μ L)	Glycogen (μ g/ml)
1	100 μ L + 0 μ L	100	2000
2	60 μ L + 40 μ L	100	1200
3	40 μ L + 60 μ L	100	800
4	20 μ L + 80 μ L	100	400
5	10 μ L + 90 μ L	100	200
6	5 μ L + 95 μ L	100	100
7	1 μ L + 99 μ L	100	20
8	0 μ L + 100 μ L	100	0

Notice : Transfer 10 μ L standard and samples into separate wells of a clear flat-bottom microplate. If the sample contains glucose, transfer an additional 10 μ L sample to another well for the Sample Blank.

3. Transfer 20 μ L Reagent C to each Sample, dH₂O Blank and Standard well. Tap plate to mix. If the sample contains glucose, transfer 20 μ L dH₂O to the Sample Blank well instead of Reagent C. Incubate 5 min at room temperature.
4. Working Reagent. For each reaction well, mix 50 μ L Reagent

产品仅供科学研究 禁止用于临床诊断、治疗
For Research Use Only. Not For Use In Diagnostic
Procedures.

A, 50 μ L Reagent B, and 2 μ L Reagent D. Transfer 100 μ L Working Reagent into each Standard and Sample well (including Sample Blank well). Tap plate to mix.

4. Incubate 5-10 min at room temperature. Read optical density at 570 nm (550-585 nm).

Fluorimetric Procedure:

For fluorimetric assays, the linear detection range is 0.4 to 200 μ g/mL glycogen. Follow steps of the colorimetric procedure, but prepare 0, 5, 10, 15 and 20 μ g/mL Standard and use a black flat-bottom microplate. Incubate 5-10 min at room temperature and read fluorescence at λ_{ex} = 530nm and λ_{em} = 585 nm.

Calculation

Subtract Blank reading (OD570nm or fluorescence intensity) from the standard reading values and plot the DOD or DF against standard concentrations. Determine the slope and calculate the glycogen concentration of the sample.

$$\text{Glycogen} = \frac{R_{\text{SAMPLE}} - R_{\text{BLANK}}}{\text{Slope}} \mu\text{g/mL}$$

R_{SAMPLE} and R_{BLANK} are the OD570nm or fluorescence intensity values of the sample and blank (water, or sample blank).

GENERAL CONSIDERATIONS

1. This assay is based on a kinetic reaction, the use of a multi-channel pipettor for adding the working reagent is recommended.
2. SH-group containing reagents (e.g., DTT, β -mercaptoethanol) may interfere with this assay and should be avoided in sample preparation.

Storage conditions

The kit is shipped on dry ice. Store all components at -20°C.

Shelf life

6 months after receipt.

Glycogen Standard Curves for Reference

