

Goat Anti- Rabbit IgG (H+L)-PE Conjugated

Technical literature is available at: www.mesgenbio.com. E-mail MesGen Technical Services if you have questions on use of this system: tech@mesgenbio.com

MESGEN
INNOVATION BIOTECHNOLOGY

Catalog Number : MAN4524

Lot Number : Refer to vial

Packaging Size : 100µL

Background

The antibody reacts with whole molecule Rabbit IgG. No antibody was detected against non-immunoglobulin serum proteins.

Conjugation

FITC. Ex: 488nm, Em: 575nm

Overview

Description : Goat polyclonal Secondary Antibody to Rabbit IgG - H&L (PE)

Reactivity : Rabbit

Formulation : Liquid

Tested applications : FC (Flow Cytometry), IF (Immunofluorescent)

IF: 1/50 ~1/1000

Since applications vary, the appropriate dilutions must be used at an assay dependent concentration.

Properties

Immunogen : Rabbit IgG whole molecule

Clonality : Goat Polyclonal

Isotype : Goat Polyclonal IgG

Storage instruction

Store at +4°C. and protected from prolonged exposure to light. Do not freeze.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Distributor in China (Mainland)

Company : Shanghai hongsheng Biotechnology Co.,Ltd

Address : Room 2-205, Xingzhan building, No.185, Liuying Road, Zhabei District, Shanghai

Tel : 86-21-56620378 **Fax** : 86-21-56620378

Email : sales@hongshengbio.com

Website : www.hongshengbio.com