

Recombinant Mouse Epidermal Growth Factor/EGF (C-6His)

MESGEN
INNOVATION BIOTECHNOLOGY

Cat.No # MGF9825

Lot.No : Refer to Vial



Packaging: 50µg \ 500µg \ 1mg

Description

EGF is a single-pass type I membrane protein, containing 8 LDL-receptor class B repeats and 9 EGF-like domains. EGF results in cellular proliferation, differentiation, and survival. EGF is a low-molecular-weight polypeptide first purified from the mouse submandibular gland, but since then found in many human tissues including submandibular gland, parotid gland. Salivary EGF, which seems also regulated by dietary inorganic iodine, also plays an important physiological role in the maintenance of oro-esophageal and gastric tissue integrity. The biological effects of salivary EGF include healing of oral and gastroesophageal ulcers, inhibition of gastric acid secretion, stimulation of DNA synthesis as well as mucosal protection from intraluminal injurious factors such as gastric acid, bile acids, pepsin, and trypsin and to physical, chemical and bacterial agents.

Amino acid sequence / Origin

Recombinant Mouse epidermal growth factor is produced by our E.coli expression system and the target gene encoding Asn977-Arg1029 is expressed with a 6His tag at the C-terminus.

Quality control

Greater than 95% as determined by reducing SDS-PAGE.

Dissolution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Store condition

Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks.

Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Aliquots of reconstituted samples are stable at < -20°C for 3 months.

For Research Use Only. Not for use in diagnostic procedures.

